

IMAN AFSHAR




SYNGO •

• SAFETY


• PHYSICS

TECHNIQUE •



- SAFETY





SAFETY

01

PATIENT SAFETY

- *Hearing*
- *SAR*
- *Clothing*
- *Frist level*
- *Laser*
- *Screening of patient*
- *Accessories*

02

MAGNET FIELD

- *5 Gauss Line*
- *ERDU*
- *RF Transmit mode*
- *GRADIENT*
- *RF*

03

STOP BOTTOM

- *TSB*
- *EPO*
- *MAGNET OFF*

04

EMREGENCY

- *Donot...*
- *Donat...*

SAFETY

01

PATIENT SAFETY

- *Hearing*
- *SAR*
- *Clothing*
- *Frist level*
- *Laser*
- *Screening of patient*
- *Accessories*

02

MAGNET FIELD

- *5 Gauss Line*
- *ERDU*
- *RF Transmit mode*
- *GRADIENT*
- *RF*

03

STOP BOTTOM

- *TSB*
- *EPO*
- *MAGNET OFF*

04

EMREGENCY

- *Donot...*
- *Donat...*

SAFETY

01

PATIENT SAFETY

- *Hearing*
- *SAR*
- *Clothing*
- *Frist level*
- *Laser*
- *Screening of patient*
- *Accessories*

02

MAGNET FIELD

- *5 Gauss Line*
- *ERDU*
- *RF Transmit mode*
- *GRADIENT*
- *RF*

03

STOP BOTTOM

- *TSB*
- *EPO*
- *MAGNET OFF*

04

EMREGENCY

- *Donot...*
- *Donat...*

SAFETY

01

PATIENT SAFETY

- *Hearing*
- *SAR*
- *Clothing*
- *Frist level*
- *Laser*
- *Screening of patient*
- *Accessories*

02

MAGNET FIELD

- *5 Gauss Line*
- *ERDU*
- *RF Transmit mode*
- *GRADIENT*
- *RF*

03

STOP BOTTOM

- *TSB*
- *EPO*
- *MAGNET OFF*



04


EMREGENCY


- *Donot...*
- *Donot...*
- *Donot...*


5gauss line


- *2.5m XY*
- *4m Z*













Spatial gradient of the static magnetic field B_0




Product of the static magnetic field B_0 and the spatial gradient of B_0




ERDU

ERDU




ERDU




EMERGENCY RAMP DOWN UNIT

QUENCH

Only to be used ...




- (1) Magnet Stop switch
- (2) Emergency Shutdown switch
- (3) Table Stop button



- (1) Magnet Stop switch
- (2) Emergency Shutdown switch
- (3) Table Stop button



Magnet Stop switch



Emergency Shutdown




Table Stop on the patient
table





Table Stop/Shutdown on
the intercom



- ***SAR IS ...***



- ***SAR IS ...***



• SAR IS ...



SARLimit(s) Exceeded

Step 1 - se_15b130

BORE Temperature 69.8 °F

Operating Mode

Current Measurement First Level

Next Measurement ☐ Normal Mode ☒ First Level

Protocol Parameters

	Current	Recommended
☒ Reduce Slices	20	16
☒ Increase TR	1340 ms	1630 ms
☒ Reduce Flip Angle	130 deg	deg

OK Open Protocol Skip

- SAR IS ...**

Patient Registration

Program Selection

Emergency Patient

Register Patient

Patient Query

Delete Procedure

* Last Name

First Name

Title

Suffix

Middle Name

* Patient ID

* Date of Birth

Mid/yyyy

Mid/yyyy

* Age

Year(s)

* Sex

* Height

ft

in

* Weight

lbs

U.S.

Admitting Diagnosis

Alerts

Medical Alerts

Allergies

Institution

Institution Name

Performing Physician

Referring Physician

Requesting Physician

Operator

Examination Information

Procedure

Accession Nr.

Reg. Proc. ID

Study Description

Study Comment

Program Selection

Select

Load Program to Queue

RF Transmit Mode

Any Polarization

Any Polarization

P

Circular Polarization only (CP)

Select

Defines the RF transmit mode for this examination.

Any Polarization means that measurement can be performed with any polarization (elliptic, circular, and linear polarization).

Circular Polarization means that measurement steps will be performed with circular polarization only (acc. to IEC 60601-2-33 ed. 3.2).

Safety relevant information needs to be validated and confirmed

* Mandatory information

Save

Cancel

Delete

Local Data

Prior Studies

Activate Windows
Go to Settings to activate Windows.

- ***SAR IS ...***

HIGH/WEIGHT

- ***SAR IS ...***

LOOKAHEAD MONITOR

HIGH/WEIGHT

- ***SAR IS ...***

RUN TIME MONITOR

LOOKAHEAD MONITOR

HIGH/WEIGHT

- ***SAR IS ...***

RUN TIME MONITOR

LOOKAHEAD MONITOR

HIGH/WEIGHT

TR/FA/#SLICE/ETL

- ***SAR IS ...***

RUN TIME MONITOR

LOOKAHEAD MONITOR

HIGH/WEIGHT

FEET FRIST

TR/FA/#SLICE/ETL

- ***SAR IS ...***

RUN TIME MONITOR

LOOKAHEAD MONITOR

HIGH/WEIGHT

FRIST LEVEL



FEET FRIST

TR/FA/#SLICE/ETL

- *SED IS ...*

- *SED IS ...*

SAR*TIME


- *SED IS ...*

SAR*TIME

An orange, seven-pointed star shape with a thin black outline, containing the text 'SAR= SPEED' in white.

SAR= SPEED

- *SED IS ...*

Two overlapping star-shaped callouts, one light orange and one medium orange, positioned below the main text. The left star contains the text 'SAR= SPEED' and the right star contains 'SED= MEATUSE'.

SAR= SPEED

SED= MEATUSE

- *SED IS ...*

SAR= SPEED

SED= MEATUSE

SAR
SED



- *SED IS ...*

SAR= SPEED


SED= MEATUSE

SAR ↓
SED ↑

14400 J/KG
6000 J/KG


General MR safety measures

- The subject must not be positioned as shown left.
- Use of paper/linen to cover all exposed skin.
- No synthetic blankets.
- No direct patient contact to magnet bore or cables.



General MR safety measures

- The subject must not be positioned as shown left.
- Use of paper/linen to cover all exposed skin.
- No synthetic blankets.
- No direct patient contact to magnet bore or cables.



- *Terminology for accessories*

01

02

03

- *Terminology for accessories*



02

03

- *Terminology for accessories*

MR Unsafe



MR Conditional



03

- *Terminology for accessories*

MR Unsafe



MR Conditional



MR

MR Safe





Implants susceptible to
electromagnetic effects



Implants susceptible to
electromagnetic effects



Open flames, no smok-
ing



Implants susceptible to
electromagnetic effects



Metallic implants and
other metallic objects
inside the body



Open flames, no smok-
ing



Implants susceptible to
electromagnetic effects



Metallic implants and
other metallic objects
inside the body



Fire extinguishers with
magnetizable metal
housing



Open flames, no smok-
ing



Implants susceptible to
electromagnetic effects



Metal parts, e.g. tools



Metallic implants and
other metallic objects
inside the body



Fire extinguishers with
magnetizable metal
housing



Open flames, no smok-
ing



Implants susceptible to electromagnetic effects



Metal parts, e.g. tools



Metallic implants and other metallic objects inside the body



Mechanical watches and electronic data carriers



Fire extinguishers with magnetizable metal housing



Open flames, no smoking



Implants susceptible to
electromagnetic effects



Metal parts, e.g. tools



Metallic implants and
other metallic objects
inside the body



Mechanical watches
and electronic data car-
riers



Fire extinguishers with
magnetizable metal
housing



Electronic data carriers



Open flames, no smok-
ing

- *LASER*



- *GRADIENT*





MRI ACCESS ZONES

ZONE

1

ZONE

2

ZONE

3

ZONE

4

ZONE
1

ZONE
2

ZONE
3

ZONE
4

ZONE
1

ZONE
2

ZONE
3

ZONE
4

ZONE

1

ZONE

2

ZONE

3

ZONE

4

ZONE

1

ZONE

2

ZONE

3

ZONE

4

STAFF LEVELS

MRI ALWAYS ON

MRI ZONES

BIOLOGICAL EFFECT

CHECKING FOR FERROMAGNETIC MATERIAL

DONOT GUESS

ANNUAL SCREENING


PREGNANT EMPLOYEES

MRI SAFETY OFFICER

FOUR IMPORTANT FIELD IN MRI


CLAUSTROPHOBIA

CONTRAST MEDIA



1

Patient Screening: Conduct thorough screening to identify any implants or conditions that may pose risks.




1

Patient Screening: Conduct thorough screening to identify any implants or conditions that may pose risks.

2

Training: Ensure all staff receive comprehensive MRI safety training.



1


Patient Screening: Conduct thorough screening to identify any implants or conditions that may pose risks.

2


Training: Ensure all staff receive comprehensive MRI safety training.

3

Magnetic Field Awareness: Be mindful of the strong magnetic field to prevent accidents with ferromagnetic objects




SAR Management: Monitor and manage Specific Absorption Rate (SAR) to minimize heating risks



SAR Management: Monitor and manage Specific Absorption Rate (SAR) to minimize heating risks



Emergency Preparedness: Have a clear protocol for emergencies like quenching or patient distress.



SAR Management: Monitor and manage Specific Absorption Rate (SAR) to minimize heating risks



Emergency Preparedness: Have a clear protocol for emergencies like quenching or patient distress.



Patient Comfort: Use tools and communication to reduce patient anxiety and improve their experience.

Safety Information

Daily MR components functionality checklist

Before using the MR system, it's always important to observe the trained safety measures and check the functionality and/or cleanliness of the following parts and areas:

- ☐ Check that all warning symbols and signs are present inside and outside the examination room.
- ☐ Check the examination, control, and equipment room for liquid spills and puddles on the floor.
- ☐ Ensure that **no magnetizable materials or objects** are present in the examination room.
- ☐ Visually inspect the table appearance and cleanliness.
- ☐ Check the Intercom functionalities (e.g., LED, sound, patient communication).
- ☐ Verify that pressing the squeeze bulb triggers the patient alert.
- ☐ Test the Table Stop buttons to ensure they function properly and stop the table immediately (see page 6).
- ☐ Check that the patient video monitoring systems transmit images properly.
- ☐ Look for the System Check icon on the Home Screen. If it's red, review the **Message History - MR Scanner and Acquisition** and contact Siemens Healthineers Service or wait until the system is ready.



Activate Win
Go to Settings t

Safety Information

Pre-screening of patients and workers

Before any MR examination, **screen and inform the patient about the effects and risks** (e.g., stimulation, acoustic noise, heating effects, patient table movement). Additionally, educate the patient on how to collaborate. Document this information exchange in a Patient consent form, which must be signed before the examination.



Ensure that patients **always remove** ferromagnetic objects before entering the examination room.

Ferromagnetic objects are:

- Bras (with clasps, underwire)
- Belts
- Trousers with zippers or buttons
- Hair pins, barrettes, or hair clips
- Removable dental braces
- Nose and mouth masks with metal strips



Safety Information

Patient safety checklist



Always prioritize patient and staff safety while performing an examination.

Before starting scanning, please check the following:

- ☐ Does the patient have any implants, are they pregnant or are they wearing any (permanent) make-up?
- ☐ Is there a possibility of embedded metal splinters (e.g., due to military activities)?
- ☐ Is the patient wearing sufficient hearing protection?
- ☐ Is the patient correctly positioned to avoid skin-to-skin contact between body parts (e.g., touching hands, feet, etc.)?
- ☐ Is the patient correctly positioned to avoid contact between the body parts and the scanner bore?
- ☐ Does the patient need active and extra medical supervision during the exam?


This is required for the following patient groups:

- Patients not capable of communicating autonomously (infants, children, patients who are confused, sedated, or unconscious)
- Decompensated cardiac patients, febrile patients, or all patients with potential of developing cardiac arrest, seizures, or claustrophobic reactions.


MR Coil Handling, Storage and Cleaning

Tips and Guidelines

Version 2022.04




Template Effective Date: 22 Nov 2021 | HILS 2223
Effective Date: 14 Jul 2022 | HOOD05162003286712
Unrestricted © Siemens Healthineers, 2022



1

Damaged RF coils are the cause of many costly and time-consuming system issues.




1

Damaged RF coils are the cause of many costly and time-consuming system issues.

2

Coil connectors/Levers/Coil plug or coil link damage



1

Damaged RF coils are the cause of many costly and time-consuming system issues.



2

Coil connectors/Levers/Coil plug or coil link damage


3


Check the coil connectors and links weekly
















Cleaning Recommendations

Cleaning Recommendations

Cleaning Recommendations

- ✓ *Clean at least every two months*
- ✓ *Prior to cleaning, switch off the LCD monitor/video display and disconnect the main power plug. Wait until it has cooled down before cleaning*
- ✓ *If additional cleaning is needed, spray window cleaner lightly on a cloth to clean the monitor. Never spray window cleaner directly on the monitor*
- ✓ *Avoid scratching the surface of the LCD monitor/video display.*

Cleaning Recommendations

Cleaning Recommendations



Approved Detergents	Not Approved Detergents
Aldehydes	Alcohol solution (70% isopropyl alcohol and 30% water) – Pads and accessories
Quaternary compounds	Benzine undiluted
Guanidine derivatives	Acetone
Peroxide compounds	Phenols
Pyridine derivatives	
Chloro derivatives	
Alcohol solution (70% isopropyl alcohol and 30% water) – COILS only	
Alcohol solution min. 60% ABV (only for coils and their accessories)	
Commercially available cleaning agents, detergent substances	
Alkylamine	
Organic acids	

Cleaning Recommendations

- ✓ *Use commercially available cleaning and disinfectant solutions. Always refer to the lists of approved detergents provided in the operator manual for each system component and accessory.*
- ✓ *Do not use hard or sharp objects (for example, knives or tweezers) to remove residue*
- ✓ *Do not pour cleaning fluid directly onto surfaces or spray with fluids. Always use a damp cloth for cleaning. Carefully wipe the plugs and connectors with a cloth. Do not touch the contacts.*
- ✓ *Do not immerse in cleaning or disinfectant liquid. Do not rinse with water*
- ✓ *Wash the straps of the patent table at a temperature of 60°C with the approved detergents*

Artifacts

Brief Overview



Template Effective Date: 27 MAY 2019 | HOOD05162003010284

Job Aid Effective Date: 19 JUL 2021 | HOOD05162003208932

Unrestricted © Siemens Healthineers 2021

Activate Win

Go to Settings to

- ✓ *In general, the image quality (IQ) depends on the image characteristics listed below. During Quality Assurance (QA), the user evaluates if the image matches the quality requirements regarding these characteristics.*
- ✓ *Artifacts are changes of the signal intensity in an MR image that do not correspond to the real tissue signal. They result mainly from physiological as well as system-related effects. There are several types of artifacts, e.g., motion artifacts, which can affect the image for different reasons.*

Characteristics	Short Description
Contrast	Differentiation between tissues with different compositions
Resolution	Image resolution is defined by the voxel size
Signal-to-Noise Ratio (SNR)	Ratio between the signal of the area of interest and the background noise
Artifact	Signal variations that can influence image appearance and interpretation

Motion Artifacts Artifacts due to Movement



- voluntary movements
- Involuntary movements

☐ ***Phase Encoding Direction***

☐ ***Ghosting artifact***

☐ ***Pulstion artifact***

- ***Advise the patient not to move during the examination.***
- ***Increase the patient's comfort if possible***
- ***repeat the scanning step.***
- ***fast or motion-insensitive mode***
- ***hyosine***






Flow Artifacts Artifacts due to Flow

- periodic repetitive motions
- involuntary movements of fluids in the body

- ❑ *Phase Encoding Direction*
- ❑ *Pulsion artifact vessel & csf*

- *change the Phase Encoding Direction*
- *FC*
- *SATURATION BAND*
- *3D*
- *SE*
- *ECG*
- *GRAPPA*



Susceptibility and Metal Implants Artifacts

- *due to magnetic sensibility relating to something*
- appear as localized image distortions or signal changes


☐ ***Generic susceptibility artifact***

☐ Metal artifacts


- higher bandwidth as well as lower TE are less affected.
- FSE VS GRE
- ETL // S.T // BASE RES // SEMAC // VAT



Generic susceptibility artifact




Metal implant artifact



Flash
BW 150 HZ/Pixel
TE 4.6 ms



Flash
BW 150 HZ/Pixel
TE 11 ms




Aliasing Artifacts Artifact due to Overfolding

- when the object to be scanned is bigger than the Field of View (FoV)

☐ *Phase Encoding Direction*

- *oversampling in Phase Direction (2D images) or Slice Direction (3D images).*
- *SATURATION BAND*
- *PHASE ENCODING DIRECTION SWAP.*




Inhomogeneity Artifacts Artifacts due to Spatial Signal Inhomogeneity

Reasons for field inhomogeneity artifacts are:


- **The radio-frequency (RF) field of the transmitting coil is not homogeneous. As a result, the RF pulses are of different strengths and the excitation is different.**
- **The radio-frequency (RF) field of the receiving coil is not homogeneous (different sensitivity). As a result, the signal is lower in less sensitive areas and stronger in areas with higher sensitivity.**
- **Inhomogeneity of the B0 magnetic field, e.g., due to shimming, eddy currents, local B0 inhomogeneities**

➤ **CALL TO SIEMENS**



Hardware Artifacts & ELECTROMAGNETIC ARTIFACT ***Artifacts due to Coil Positioning, Broken or Wrong Coil***

- ***First, check that the coil is positioned correctly.***
- ***A broken coil element can also lead to a signal loss and increase in noise.***
- ***SHAPE & DISTANCE & SIZE***



ساعت:

دمای اتاق ها			تست Shimming با فانوم	سیستم چیلر			گیج های کابینت Cooling			مقدار Helium	Alarm Box عملکرد اتاق کنترل (سبز بودن تمام چراغها)	
کنترل	نگینگر	مگنت	Ok Not Ok	دبی آب ورودی چیلر به کابین Lit/Min (110-120) 11-12 rpm	فشار آب ورودی چیلر به کابین (bar) (1.5 - 5)	دمای آب ورودی چیلر به کابین (سانتigrاد) (8-13)	فشار گیج A (bar) 3< بالاتر از ۳	کارکرد کمپرسور (ساعت) کمتر از ۳۰۰۰۰ هر 30 هزار ساعت فیلتر عوض شود)	فشار گاز کمپرسور (bar) (19-22)		کنترل روزانه حجم هلیوم با ذکر عدد از قسمت cooling status از روی پنل ۱۵۰۰ قابل قبول است تحت نام پائول بیشتر از این عدد گزارش شود.	Ok

توضیحات	تجهیزات ام آر آی	عملکرد کویل ها
	اکسیرن دیواری	
	بالس اکسیمتر (نحوه استقرار و شارژ مناسب)	
	عملکرد سرور سیستم پیام رسان هوشمند	
	خطاهای روزانه Event Log Tab دستگاه	
	کنترل فضای Database (بالای ۶۰ درصد پر نباشد)	
	فعال بودن Auto Send تصاویر به پکس و ارسال تصاویر روز قبل	
	PMU (Physiologic Monitoring Unit)	
	Breast	
	Large Loop	
	Small Loop	
	Large Flex	
	Small Flex	
	Body 2	
	Body 1	
	Knee	
	Shoulder	
	Cervical	
	Brain	

وضعیت		نام قسمت ها
☐ مناسب	☐ نامناسب	عملکرد گیت ورودی اتاق مگنت
☐ سالم	☐ ناسالم	عملکرد مانیتور ام آر ای در اتاق کنترل
☐ سالم	☐ ناسالم	عملکرد گوشی بیمار در اتاق مگنت
☐ وجود دارد	☐ وجود ندارد	وجود ویلچر و تخت مخصوص ام آر ای
☐ سالم	☐ ناسالم	عملکرد اینتر کام (ارتباط کلامی با بیمار)
☐ مناسب	☐ نامناسب	بنجره ها و صندلی ها
☐ مناسب	☐ نامناسب	کف پوش ها و سالنها
☐ مناسب	☐ نامناسب	عدم وجود خون یا زائده فلزی در سوکت کوبلهای روی تخت
☐ مناسب	☐ نامناسب	رختکن ها
☐ مناسب	☐ نامناسب	ریکاوری، تختها و تخت ام آر ای
☐ مناسب	☐ نامناسب	روشویی و تی شویی و کنترل طی ها طبق رنگ بندی
☐ مناسب	☐ نامناسب	اتاقهای تکنیک، مگنت، کنترل
☐ انجام شده	☐ انجام نشده	تفکیک پسماندهای بیمارستانی
☐ وجود دارد	☐ وجود ندارد	وجود لوازم مصرفی پزشکی و اداری کافی (CD/DVD - لوازم حفاظت فردی - لوازم حفاظت از بیمار دروشیت، گان)
☐ وجود دارد	☐ وجود ندارد	وجود مواد ضد عفونی کننده در ریکاوری، اتاق کنترل و مگنت
☐ مناسب	☐ نامناسب	(کنترل تاریخ و حجم سطل) Safety Box
☐ انجام شد	☐ انجام نشده	تکمیل چک لیست ترالی کد و ملزومات متعلقه (اکسیژن - داروها - دفیلیراتور - پالس - مانیتور - ساکشن)
توضیحات/ اقدامات:		