


Head of Ventilator






Head Connections to Patients

1. Connection for the expiratory set.
2. Connections for pneumotachograph hoses of the expiratory valve.
3. O2 cell
4. Connection for the nebulization hose.
5. Connection Towards Patient of the patient circuit.
6. Inlet for the capnograph connector.
7. Connection for proximal pneumotachograph.





How to remove expiratory valve





How to maintain expiratory valve

20 min
70% alcohol







O₂ sensor

Adjusted FiO₂ should be equal to monitored FiO₂

- Defectiveness
- Low O₂ Pressure
- Incorrect O₂ inlet connection



Back of Ventilator

Speaker

Fan

ON/OFF

Safety
Clamp






Compressor

- Make compressed medical air
- Filtered and Prepared to mixed with O2
- Time Service : Evry 4000 Hours
- Pressure range: **3.5 – 7 bar**



Back of Compressor



Change every 6
months



Humidifier



- **Non-Servo (For ADL)**





Humidifier

- **Non-Servo (For ADL)**
 - **How to Connect**



Just fill with Distilled Water
Do not use with HME Filter simultaneously



Humidifier

- **Non-Servo (For ADL)**
 - Chamber Maintenance
- **Autoclavable**
 - 120 °C (393 °F)
 - 96 kPa
 - 15 min
- **Eto**
 - 55 °C (131 °F)





HME filter (Heat and moisture exchanger)

- **Contains alkaline salts**
- **Just warm and humidify the inspiration**
- **Best installation : in front of Y piece**
- **Do not confuse this with antibacterial filter**
- **Make sure the filter is not blocked**
- **Installation in expiratory valve:**
 - Returning excess air to the lung**
- **If you sense any change in filter color or resistance in airway, you must change the filter**



Overall View





How to operate



Initial Test

- Nebulizer should not connect to device




Initial Test

- **Active Humidifier**
Press Knob to start




Initial Test

- **Passive**
Press Knob to start



Initial Test

- Calibrating
O₂ sensor
Expiratory Valve
- Check leaks
- Calculating dead space and airway resistance



Initial Test Done choose your mode



ADL

| L/min Peak Flow | s Ii | s Ie | Current Ib | bpm f _{total} | L V _I | L/min V̇ _E | % Oxygen |
|--------------------|---------|---------|---------------|---------------------------|---------------------|--------------------------|-------------|
| 30 | 0.52 | 0.32 | 12 | 3 | 60 | 40 | |

Alarm Signals

Vent Inop

High

Med Low

Alarm Limits

30-60 s

Menu

Mode

Volume compensation



Ventilation Modes

Mode PCV

Press control knob to start

Press Esc to exit

| | | | | | | | |
|--------------------|--------------|------------|----------------|--------------------|--------------------|-------------------|--------------|
| 15 cmH2O PCV | 1 e Ti | 1:4 I:E | 12 bpm f | 5 cmH2O PEEP | 3 L/min Sens | 50 % Oxygen | Rise Time |
|--------------------|--------------|------------|----------------|--------------------|--------------------|-------------------|--------------|

110-220 VAC

Battery/charge



Suction %O₂

Manual Insp.

Insp/Exp Pause

Ctrl

Esc

EDP-TS

P2 P1 Prox

O₂%

CO₂

ADL

PCV

| L/min Peak Flow | s Ii | s Ie | Current Ib | bpm f _{total} | L V _I | L/min V̇ _E | % Oxygen |
|--------------------|---------|---------|---------------|---------------------------|---------------------|--------------------------|-------------|
| 75.1 | 1 | 3.22 | 1:3.1 | 14 | 0.564 | 7.9 | 50 |

Alarm Signals

Vent Inop

High

Med Low


Alarm Limits

30-60 s


Menu

Mode

Volume compensation



High tidal volume



| | | | | | | | |
|--------------------|--------------|------------|----------------|--------------------|--------------------|-------------------|--------------|
| 15 cmH2O PCV | 1 e Ti | 1:4 I:E | 12 bpm f | 5 cmH2O PEEP | 3 L/min Sens | 50 % Oxygen | Rise Time |
|--------------------|--------------|------------|----------------|--------------------|--------------------|-------------------|--------------|

110-220 VAC

Battery/charge



Suction %O₂

Manual Insp.

Insp/Exp Pause

Ctrl


Esc

EDP-TS

P2 P1 Prox

O₂%

CO₂



Indication of Use :

Ventilation

Oxygenation

Airway Protection



VENTILATION (CO₂ WASHING)

$$MV = VT \times \text{RATE}$$

- VT
- RATE
- I/E (TI OR PEAK FLOW)

Tidal Volume Estimation

- IBW

Male: $50 + 0.91(\text{Patient Height} - 152.4)$ or Patient Height – 105

Female: $45.5 + 0.91(\text{Patient Height} - 152.4)$ or Patient Height – 107

$$\text{Tidal Volume} = \text{IBW} \times (6 - 8 \text{ ml/Kg}) \quad (\text{Peak Pressure} < 30 - 40)$$

PC VC

Rate Readjustment

Initial rate (12 – 16) → set on mean (14)

ABG Test → Amount of CO₂ if > 45

$$\text{New Rate} = \frac{\text{CO}_2 \times \text{Previous Rate}}{40}$$




OXYGENATION

- FIO₂
- PEEP
- MAP

- FiO₂ 100%(40-60)% SPO₂(92-94) OR PaO₂(65-70)
- PEEP 5 cmH₂O (Physiological 3 – 5 cmH₂O)

SYNCHRONY

- TRIGGER (2-3)
- PSV (8-10)



PULMONARY DISEASES

| Type of Disease | Name of Disease | Common Symptoms |
|-----------------|--|---|
| Obstructive | COPD (Chronic Obstructive Pulmonary Disease) | Shortness of breath, chronic cough, sputum production, wheezing |
| | Asthma | Episodic shortness of breath, wheezing, cough (worse at night or early morning), chest tightness |
| | Bronchiectasis | Chronic productive cough, shortness of breath, recurrent lung infections |
| | Tracheobronchial Fistula | Chronic cough, recurrent pneumonia, abnormal breath sounds |
| Restrictive | Pulmonary Fibrosis | Progressive shortness of breath, dry cough, fatigue, cyanosis |
| | ARDS (Acute Respiratory Distress Syndrome) | Sudden severe shortness of breath, rapid breathing, low oxygen levels, often requiring mechanical ventilation |
| | Severe Scoliosis (affecting lung function) | Reduced lung capacity, exertional shortness of breath |
| | Neuromuscular Diseases (ALS, Muscular Dystrophy) | Respiratory weakness, reduced vital capacity, shallow breathing |
| | Sarcoidosis | Shortness of breath, dry cough, fatigue, systemic involvement |
| | Chronic Pleural Effusion or Pleural Fibrosis | Chest pain, shortness of breath, reduced breath sounds |

THANKYOU